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- ATLAS, R.M., L.C. PARK (1993) Handbook of Microbiological Media for the examination of Food. CRC Press Inc. Boca Ratón.
- DOWNES, F.P. & K. ITO (2001) Compendium of Methods for the Microbiological Examination of Foods. 4th ed. APHA. Washington DC. USA.
- EUROPEAN PHARMACOPOEIA 8.0 (2014) 8th ed. § 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms. Harmonised Method. EDQM. Council of Europe. Strasbourg.
- HORWITZ, W. (2000). Official Methods of Analysis of the AOAC Internacional. 17th ed. Gaithersburg Md. USA.
- ICMSF (1978) Microorganisms in Foods 1. University of Toronto Press.
- ISO 6340:1995 STANDARD. Water Quality - Detection of *Salmonella* spp.
- ISO 11133:2014/ Adm 1:2018. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- PASCUAL ANDERSON, M^a R. (1992) Microbiología Alimentaria. Díaz de Santos, S.A. Madrid.
- TAYLOR, W.J. (1965) Isolation of *Shigella*. I. Xylose Lysine Agars: New media for isolation of enteric pathogens. Am. J. Clin. Path 44:471-475.
- US FDA (Food and Drug Administrations). (1998) Bacteriological Analytical Manual. 8th ed. Revision A. AOAC International. Gaithersburg, Md. USA.
- USP 33 - NF 28 (2011) <62> Microbiological examination of non-sterile products: Test for specified microorganisms. Harmonised Method. USP Corp. Inc. Rockville. MD. USA.