

Gentamicin Cat. 6805

For E.coli in molecular genetics studies

Practical information

Aplications	Categories	
Selection of transformants	Escherichia coli	
Preparation and recovery of competent cells	Escherichia coli	

Industry: Culture media ingredients / Culture media for Molecular biology

Principles and uses

Gentamicin is used with LB Media for the isolation of bacteria modified with a plasmid that includes a chloramphenicol resistance gene.

Gentamycin inhibits protein synthesis by binding to L6 protein of the 50S ribosomal subunit.

Formula per vial

Gentamicin (g)		5	

Typical formula g/L * Adjusted and/or supplemented as required to meet performance criteria.

Preparation

Add gentamicin powder or a sterile gentamicin solution to LB medium (Cat.1551, Cat.1552, Cat.1266, Cat.1308, Cat.1231, Cat.1083) to get the desired final concentration (µg/mL). Autoclaved and cooled to 50°C. Mix well and distribute into sterile containers.

Instructions for use

The Stock solution can be prepared with water, in a concentration of 10-50 mg/ml and should be stored at -20°C.

The working concentrations is 5 µg/ml for stringent plasmids and a concentration of 20 µg/ml for relaxed plasmids.

Quality control

Solubility	Appareance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°C)
N/A	Powder	White or almost white	N/A	3.5-5.5

Physical-chemical characteristics

Description	Specification
Loss on drying	<18%
рН	3,5-5,5
Identification	Positive
Specifical optical rotation	+107° to 121°
Residue after ignition	<1%
Limit of Methanol	<1%
Content of gentamicins:	C1+C2b : 25% to 50% / C1a : 10% to 35% / C2+C2a: 25% to 55%
Assay (anhydrous basis)	>590 µg/mg

Storage

Temp. Min.:0 °C Temp. Max.:5 °C

Bibliography

Ryan N. (1985) Personal communication.
Rogol M., Sechter I., Grinberg L., Gerichter Ch. B. (1992) J. Med. Microbiol. 12. 229-231.
Atlas, R.M., L.C.Parks (1993) Handbook of Microbiological Media. CRC Press, Inc. London
The condensed protocols from molecular cloning: a laboratory manual/ Joseph Sambrook, David W .Russell