

Specification

General purpose solid medium containing animal and plant peptone according to Pharmacopoeial Harmonised Method and ISO Standards.

Presentation

20 Plates /Irradiated
90 mm - Double wrapping
with: 21 ± 2 ml

Packaging Details

1 box with 2 cellophane bags (double wrapping) with 10 plates/bag. Side labeling. Every pack exhibits a irradiation indicator stacked on the side of the bag.(8-14kGy).

Shelf Life

3,5 months

Storage

2-14 °C

Composition

Composition (g/l):

Peptone from casein	15.00
Soya peptone.....	5.00
Sodium chloride.....	5.00
Agar.....	15.00

Description /Technique

Description

TSA is a widely used medium containing two peptones which support the growth of a wide variety of organisms, even that of very fastidious ones such as Neisseria, Listeria, Brucella, etc. It is frequently used for routine diagnostic purposes due to its reliability and its easily reproducible results.

Technique

This medium can be inoculated directly or after enrichment broth.

Spread the plates by streaking methodology or by spiral method.

The inoculated plates are incubated at 30-35 °C for 24-72 h (bacteria) and 3-5 days for fungi (yeast & molds). Examined daily

(Incubation times greater than those mentioned above or different incubation temperatures may be required depending on the sample, on the specifications).

Each laboratory must evaluate the results according to their specifications.

Quality control

Physical/Chemical control

Color : Straw-coloured yellow pH: 7.3 ± 0.2 at 25°C

Microbiological control

Growth Promotion Test 50-100 CFU according to harmonized pharmacopoeial monographs and test methods & ISO 11133:2014/A1:2018

Spiral Spreading: Practical range 50 - 100 CFU (productivity).

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020.

Aerobiosis. Incubation at 30-35 °C. Read after 18-24 h to 72 h for bacteria and 3-5 days for fungi.

Microorganism

Escherichia coli ATCC® 8739, WDCM 00012

Staphylococcus aureus ATCC® 6538, WDCM 00032

Bacillus subtilis ATCC® 6633, WDCM 00003

Candida albicans ATCC® 10231, WDCM 00054

Ps. aeruginosa ATCC® 9027, WDCM 00026

Aspergillus brasiliensis ATCC® 16404, WDCM 00053

Salmonella typhimurium ATCC® 14028, WDCM 00031

L. monocytogenes ATCC® 13932, WDCM 00021

Clostridium perfringens ATCC® 13124, WDCM 00007 (37°C)

Clostridium sporogenes ATCC® 19404, WDCM 00008

Sph. aureus ATCC® 25923, WDCM 00034

Growth

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Good (≥70%)

Sterility Control

Incubation 48 h at 30-35 °C and 48 h at 20-25 °C: NO GROWTH.

Check at 7 days after incubation in same conditions.

Bibliography

- ATLAS, R.M. & L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press, Inc. London.
- COLIPA (1997) Guidelines on Microbial Quality Management (MQM). Brussels.
- DOWNES, F.P. & K. ITO (2001) Compendium of Methods for the Microbiological Examination of Food, 4th ed, ASM, Washington D.C.
- EUROPEAN PHARMACOPOEIA 10.0 (2020) 10th ed. § 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms. Harmonised Method. EDQM. Council of Europe. Strasbourg.
- FDA (Food and Drug Administrations) (1998) Bacteriological Analytical Manual. 8th ed. Revision A. AOAC International. Gaithersburg. MD.
- HORWITZ, W. (2000) Official Methods of Analysis of AOAC INTERNATIONAL, 17th ed. Gaithersburg, MD. USA.
- ISO 9308-1 Standard (2000) Water Quality. Detection and enumeration of E. coli and coliform bacteria. Membrane filtration method.
- ISO 11731 Standard (2017) Water Quality. - Enumeration of Legionella.
- ISO 11133:2014/ Adm 1:2018. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- ISO 18415 Standard (2017) Cosmetics - Microbiology - Detection of specified and non-specified microorganisms.
- ISO 21149 Standard (2017) Cosmetics - Microbiology - Enumeration and detection of aerobic mesophilic bacteria.
- ISO 21150 Standard (2015) Cosmetics - Microbiology - Detection of Escherichia coli.
- ISO 22717 Standard (2015) Cosmetics - Microbiology - Detection of Pseudomonas aeruginosa.
- ISO 22718 Standard (2015) . Cosmetics - Microbiology - Detection of Staphylococcus aureus.
- ISO 22964 (2017) Microbiology of the food chain.- Horizontal method for the detection of *Cronobacter spp*
- PASCUAL ANDERSON, M[®]R[®] (1992) Microbiología Alimentaria. Díaz de Santos S.A., Madrid.
- USP 33 - NF 28 (2011) <62> Microbiological examination of non-sterile products: Test for specified microorganisms. Harmonised Method. USP Corp. Inc. Rockville. MD. USA.